



**SINGER
INSTRUMENTS**
A RESPONSIBILITY TO SCIENCE!



Automating MALDI-TOF sample preparation
protocols to facilitate bacterial identification.

Automating MALDI-TOF sample preparation protocols to facilitate bacterial identification



APPLICATION

Michael D. Nairn¹, Philip Kirk², Matthew E. Openshaw¹, Oliver Severn², Leah Ashley²

1, Shimadzu, Manchester, UK; 2, Singer Instruments, Roadwater, UK

ABSTRACT

Matrix-assisted laser desorption/ionisation-time of flight (MALDI-TOF) is a mass spectrometry (MS) that characterises various intact biomolecules. In microbiology, MALDI-TOF aids pathogen detection and microbial profiling by identifying bacteria to the strain level. However, MALDI-TOF requires rigorous sample processing methods to ensure the reproducibility of bacterial cells from MS data. In this application, we show that introducing automation into MALDI-TOF pipelines can reduce sample processing times while maintaining accurate bacterial identification. This document provides you with the blueprint to adopt automation into your MALDI-TOF pipelines with the PIXL colony picker.

Identifying bacteria at the species level underpins clinical and environmental microbiology. The introduction of MALDI-TOF has greatly increased the resolution in bacterial profiling, largely replacing more traditional qualitative staining assays for characterising bacteria. In the clinic, MALDI-TOF is helping diagnose disease¹ and identify antibiotic-resistant bacteria². Microbiologists also use MALDI-TOF to distinguish bacterial strains and analyse human and environmental microbiomes^{3,4}. MALDI-TOF's success in conducting species-level bacterial identification lies in profiling the proteins that a microorganism produces⁶. Different bacterial species have unique protein signatures, from their ribosomal proteins to membrane proteins, that MALDI-TOF can distinguish.

Although MALDI-TOF can recognise bacterial species by its proteins, generating reproducible MS data for bacterial identification requires many steps. Variables within these steps, such as the sample preparation method, the amount of bacterial colony obtained, and the age of the colony, can affect the quality of the MS data⁷. With many experimental sources of variation present in MALDI-TOF, researchers have increasingly looked to automation to minimise technical errors.

THE QUESTION: AUTOMATION AND MALDI-TOF

Integrating automation provides several benefits for researchers. Most notably, automating experimental workflows reduces human errors, increases reproducibility, and enhances researcher efficiency⁸. For MALDI-TOF MS workflows, automating the sample processing and species identification protocols could make species-level calls more reproducible and accurate. In this application, we investigated whether the colony sample preparation procedure could be automated for MALDI-TOF MS species identification.



METHODS

An automated MALDI-TOF pipeline was developed to streamline bacterial identification from lab cultures (Figure 1).

This pipeline was tested on *Escherichia coli* for method preparation and *Bacillus subtilis* (*B. subtilis*), *Cellulomonas uda* (*C. uda*), and *Pantoea agglomerans* (*P. agglomerans*) for testing the automated method. In the automated pipeline, colonies were picked from cultured samples and deposited onto a disposable MALDI MS sample plate, such as the FlexiMass-DS slides (Shimadzu), using a PIXL precision microbial colony picker (Singer Instruments). The automated protocol was compared with manual preparations where an inoculation loop was used to spot the plates with sample material. The automation protocol was then further refined by determining whether it should feature a pinning method or a smearing method (Figure 2).

To further evaluate the automated pipeline, different areas of the culture plates were also picked for each bacterial species (Figure 3).

Once the samples were prepared, the pins and smears were analysed on an iDplus Performance MALDI-TOF mass spectrometer (Shimadzu) and submitted to a SARAMIS database for identification.




Figure 1. The PIXL colony picking robot [left] and its inner components [above], the automated system used for MALDI-TOF sample processing.




Figure 2. The two automated bacterial transfer methods employed. The pin method (far left) 'dabs' the picked sample onto a spot on the target surface, while the smear method (left) forms a smear across the target surface.




Figure 3. Three areas where bacteria were picked. Green refers to the thickest part of the culture growth, or "Deep". Pink represents the "Edges" of the culture growth streak. Yellow refers to the material picked from the thickest part of the culture growth 'steak' (solid arrows, 1) then dabbed onto an area of agar not containing any culture to remove excess material prior to depositing onto the MALDI target (dashed arrows, 2).

APPLICATION

RESULTS

To begin assessing the automated sampling method for MALDI-TOF, researchers tested it on three types of colonies across the culture plates. Doing so showed that the automated smear and automated single/double deposition of matrix solution were equally accurate in identifying *E. coli* with >99.9% confidence (Figure 4). In contrast, the manual method introduced plasticisers and other contaminants that obfuscated the *m/z* peaks for identifying microbial species.

With the automated protocol developed, researchers then determined whether the automated sample collection could identify three other bacterial species with high confidence. The automated protocol allowed two of the species, *P. agglomerans* and *C. uda*, to be identified with over 99.9% confidence. While there were red IDs for the *B. subtilis* identification, this was due to a mixed ID in the database, not the automated sampling method itself. Nonetheless, most promising about the data was the MALDI-TOF MS spectra for *P. agglomerans* being nearly identical for both the PIXL automation and manual approaches (Figure 6).



APPLICATION

	acquisition time	name	sample	%	family	genus	species	ty	sl	datac	matrix
Edge (single dip matrix)	26 Apr 2023 19:17	slide2_edge_singledep_0001_2D4p		99.90	Family 1 Enterobacter	Escherichia	coli				187
	26 Apr 2023 19:17	slide2_edge_singledep_0001_2D3p		99.90	Family 1 Enterobacter	Escherichia	coli				184
	26 Apr 2023 19:17	slide2_edge_singledep_0001_2D2p		99.90	Family 1 Enterobacter	Escherichia	coli				185
	26 Apr 2023 19:17	slide2_edge_singledep_0001_2D1p		99.90	Family 1 Enterobacter	Escherichia	coli				186
	26 Apr 2023 19:17	slide2_edge_singledep_0001_2C4p		99.90	Family 1 Enterobacter	Escherichia	coli				185
	26 Apr 2023 19:17	slide2_edge_singledep_0001_2C3p		99.90	Family 1 Enterobacter	Escherichia	coli				185
Edge (double dip matrix)	26 Apr 2023 19:17	slide2_edge_doubledep_0001_2C2p		0.00							169
	26 Apr 2023 19:17	slide2_edge_doubledep_0001_2C1p		99.90	Family 1 Enterobacter	Escherichia	coli				172
	26 Apr 2023 18:45	slide2_edge_doubledep_0001_2D4p		99.90	Family 1 Enterobacter	Escherichia	coli				231
	26 Apr 2023 18:45	slide2_edge_doubledep_0001_2D3p		99.90	Family 1 Enterobacter	Escherichia	coli				200
	26 Apr 2023 18:45	slide2_edge_doubledep_0001_2D2p		99.90	Family 1 Enterobacter	Escherichia	coli				214
	26 Apr 2023 18:45	slide2_edge_doubledep_0001_2D1p		99.90	Family 1 Enterobacter	Escherichia	coli				243
Deep + blank dab (single dip matrix)	26 Apr 2023 18:45	slide2_deep_singledep_0001_2D4p		99.90	Family 1 Enterobacter	Escherichia	coli				236
	26 Apr 2023 18:45	slide2_deep_singledep_0001_2D3p		99.90	Family 1 Enterobacter	Escherichia	coli				225
	26 Apr 2023 18:45	slide2_deep_singledep_0001_2D2p		99.90	Family 1 Enterobacter	Escherichia	coli				242
	26 Apr 2023 18:45	slide2_deep_singledep_0001_2D1p		99.90	Family 1 Enterobacter	Escherichia	coli				227
	26 Apr 2023 18:57	slide2_deep_singledep_0001_2F4p		99.90	Family 1 Enterobacter	Escherichia	coli				191
	26 Apr 2023 18:57	slide2_deep_singledep_0001_2F3p		99.90	Family 1 Enterobacter	Escherichia	coli				154
Deep + blank dab (double dip matrix)	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2F2p		99.90	Family 1 Enterobacter	Escherichia	coli				204
	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2F1p		99.90	Family 1 Enterobacter	Escherichia	coli				156
	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2E4p		99.90	Family 1 Enterobacter	Escherichia	coli				188
	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2E3p		99.90	Family 1 Enterobacter	Escherichia	coli				183
	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2E2p		99.90	Family 1 Enterobacter	Escherichia	coli				192
	26 Apr 2023 18:57	slide2_deep_doubledep_0001_2E1p		99.90	Family 1 Enterobacter	Escherichia	coli				148
Deep (single dip matrix)	26 Apr 2023 18:50	slide2_deep_singledep_0001_2L4p		99.90	Family 1 Enterobacter	Escherichia	coli				197
	26 Apr 2023 18:50	slide2_deep_singledep_0001_2L3p		99.90	Family 1 Enterobacter	Escherichia	coli				180
	26 Apr 2023 18:50	slide2_deep_singledep_0001_2L2p		99.90	Family 1 Enterobacter	Escherichia	coli				205
	26 Apr 2023 18:50	slide2_deep_singledep_0001_2L1p		99.90	Family 1 Enterobacter	Escherichia	coli				208
	26 Apr 2023 18:50	slide2_deep_doubledep_0001_2K4p		99.90	Family 1 Enterobacter	Escherichia	coli				226
	26 Apr 2023 18:50	slide2_deep_doubledep_0001_2K3p		99.90	Family 1 Enterobacter	Escherichia	coli				199
Deep (double dip matrix)	26 Apr 2023 18:50	slide2_deep_doubledep_0001_2K2p		99.90	Family 1 Enterobacter	Escherichia	coli				210
	26 Apr 2023 18:50	slide2_deep_doubledep_0001_2K1p		99.90	Family 1 Enterobacter	Escherichia	coli				219
	26 Apr 2023 19:11	slide2_deep_doubledep_0001_2B4p		99.90	Family 1 Enterobacter	Escherichia	coli				176
	26 Apr 2023 19:11	slide2_deep_doubledep_0001_2B3p		99.90	Family 1 Enterobacter	Escherichia	coli				185
	26 Apr 2023 19:11	slide2_deep_doubledep_0001_2B2p		90.20	Family 1 Enterobacter	Escherichia	coli				123
	26 Apr 2023 19:11	slide2_deep_doubledep_0001_2B1p		99.90	Family 1 Enterobacter	Escherichia	coli				240
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2A4p		99.90	Family 1 Enterobacter	Escherichia	coli				183
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2A3p		99.90	Family 1 Enterobacter	Escherichia	coli				138
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2A2p		99.90	Family 1 Enterobacter	Escherichia	coli				148
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2A1p		94.00	Family 1 Enterobacter	Escherichia	coli				148
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2H4p		99.90	Family 1 Enterobacter	Escherichia	coli				207
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2H3p		99.90	Family 1 Enterobacter	Escherichia	coli				194
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2H2p		99.90	Family 1 Enterobacter	Escherichia	coli				215
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2H1p		99.90	Family 1 Enterobacter	Escherichia	coli				225
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2G4p		99.90	Family 1 Enterobacter	Escherichia	coli				204
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2G3p		99.90	Family 1 Enterobacter	Escherichia	coli				225
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2G2p		99.90	Family 1 Enterobacter	Escherichia	coli				231
	26 Apr 2023 19:05	slide2_deep_doubledep_0001_2G1p		99.90	Family 1 Enterobacter	Escherichia	coli				185

Figure 4. Identification results for *E. coli* samples prepared with the automated smear pipeline for samples taken from all areas of the culture plate. Dark green = >99.9% confidence, light green = 90.0 – 99.8% confidence, red = mixed ID, white = no ID.

27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2L4[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		222
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2L3[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		204
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2L2[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		214
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2L1[c]		95.00	Family I Enterobacteri	Pantoea	agglomerans		214
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2K4[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		218
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2K3[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		217
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2K2[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		226
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2K1[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		202
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2J4[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		175
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2J3[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		194
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2J2[c]		98.60	Family I Enterobacteri	Pantoea	agglomerans		153
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2J1[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		226
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2H4[c]		95.20	Family I Enterobacteri	Pantoea	agglomerans		180
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2H3[c]		95.20	Family I Enterobacteri	Pantoea	agglomerans		189
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2H2[c]		99.90	Family I Enterobacteri	Pantoea	agglomerans		157
27 Apr 2023 18:03	slide2_Pagg_doubledip_0001_2H1[c]		91.80	Family I Enterobacteri	Pantoea	agglomerans		147
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2H4[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		168
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2H3[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		197
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2H2[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		206
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2H1[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		187
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2G4[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		178
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2G3[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		186
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2G2[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		175
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2G1[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		178
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2F4[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		155
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2F3[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		166
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2F2[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		190
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2F1[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		189
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2E4[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		181
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2E3[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		179
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2E2[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		154
27 Apr 2023 18:19	slide2_Cuda_doubledip_0001_2E1[c]		99.90	Family VI Cellulomonad	Cellulomonas	uda		159
27 Apr 2023 18:29	slide2_Bsubt_doubledip_0001_2D4[c]		99.90	Family I Bacillaceae	Bacillus	subtilis		171
27 Apr 2023 18:29	slide2_Bsubt_doubledip_0001_2D3[c]		99.90	Family I Bacillaceae	Bacillus	subtilis		155

Figure 5. Identification results for the bacterial species *B. subtilis*, *C. uda*, and *P. agglomerans* when using the optimised automated method. Dark green = >99.9% confidence, light green = 90.0 – 99.8% confidence, red = mixed ID.




Figure 6. Representative MALDI-TOF MS spectra for *P. agglomerans* with a manually prepared smear sample (in red) and a smear sample prepared with the automated protocol (in blue).

CONCLUSIONS

Integrating automation into research pipelines can streamline research protocols and improve reproducibility in biomedical research. In microbiology, automation may also prove useful for identifying microbial species with MALDI-TOF. Here, we show that automation can facilitate sample collection and streaking protocols, both of which affect the accuracy and reproducibility of MALDI-TOF MS-based bacterial identification.

Using the PIXL to automate colony picking and matrix deposition allowed reference microbial species to be identified correctly and reproducibly. Furthermore, the automated colony-picking method generated MALDI-TOF MS spectra like those produced with the manual method. Put together, automating the sample collection protocol makes MALDI-TOF identification more reproducible.

With automated MALDI-TOF MS pipelines now possible, researchers can efficiently type and characterise microbes for a wide range of applications, from pathogen identification to antibiotic sensitivity testing.

Contact us at contact@singerinstruments.com to learn more about the [PIXL colony picker](#) and how it can be seamlessly integrated into your bacterial identification workflows.



APPLICATION

REFERENCES

1. Patel, R. (2015) MALDI-TOF MS for the Diagnosis of Infectious Diseases. *Clinical Chemistry*, 61(1), pp. 100-111. doi: 10.1373/clinchem.2014.221770.
2. Kostrzewa, M., Sparbier, K., Maier, T., & Schubert, S. (2013) MALDI-TOF MS: an upcoming tool for rapid detection of antibiotic resistance in microorganisms. *PROTEOMICS – Clinical Applications*, 7(11-12), pp. 767-778. doi: 10.1002/prca.201300042.
3. Mojumdar, A., Yoo, H. J., Kim, D. H., & Cho, K. (2022) High-Throughput Screening Technique for Microbiome using MALDI-TOF Mass Spectrometry: A Review. *Mass Spectrometry Letters*, 13(4), pp. 106-114.
4. Ashfaq, M. Y., Da'na, D. A., & Al-Ghouti, M. A. (2022) Application of MALDI-TOF MS for identification of environmental bacteria: A review. *Journal of Environmental Management*, 305, 114359. doi: 10.1016/j.jenvman.2021.114359.



SINGER INSTRUMENTS

A RESPONSIBILITY TO SCIENCE!

Roadwater
Watchet
Somerset
TA23 0RE, UK.

+44 (0)1984 640226 (tel)

+44 (0)1984 641166 (fax)

contact@singerinstruments.com

singerinstruments.com



**Scan to visit the
PIXL webpage.**

OUR DISTRIBUTORS

For a list of our distributors visit:

singerinstruments.com/friends



/SingerInstruments